

Masterflex® Platinum-Cured Silicone Pump Tubing


Resist heat



Applications: Pharmaceutical, biotech, peristaltic pumping, and food applications.
Characteristics: Excellent heat resistance, low compression set, odorless, nontoxic. Translucent.
Certifications: USP Class V (Exceeds Class VI, European Pharmacopoeia (EP 3.1.9), FDA 21 CFR 177.2600, FDA 21 CFR 210 and 211, exceeds 3A, cGMPs, REACH (non-DEHP) Compliant, RoHS Compliant, ADFC Compliant.
Temperature range: -58 to 446°F (-50 to 230°C)
Sterilize: Autoclave or gamma irradiation.

Masterflex® BioPharm Plus Silicone Pump Tubing

Get longer life



Applications: Pharmaceutical, cosmetic, research, and autoanalysis applications requiring long service life and lower gas permeability.
Characteristics: Platinum-cured silicone. Ultra-smooth inner liner with extremely low leachables. Lower spallation and excellent pumping service life. Translucent.
Certifications: USP Class VI, European Pharmacopoeia (EP 3.1.9), FDA 21 CFR 177.2600, exceeds 3A, REACH (non-DEHP) Compliant, RoHS Compliant.
Temperature range: -75 to 450°F (-60 to 232°C)
Sterilize: Autoclave, ethylene oxide, or gamma irradiation.



ANIMAL DERIVED
INGREDIENT FREE

Tubing size	Hose barb in. (mm)	Pack			Spool		
		Catalog number	Ft/ pk	Price/ pk	Catalog number	Ft/ spool	Price/ spool
L/S® Precision pump tubing							
L/S 13	1/16 (1.6)	GY-96410-13	25		GY-96403-13	500	
L/S 14	1/16 (1.6)	GY-96410-14	25		GY-96403-14	500	
L/S 16	1/8 (3.2)	GY-96410-16	25		GY-96403-16	500	
L/S 25	3/16 (4.8)	GY-96410-25	25		GY-96403-25	500	
L/S 17	1/4 (6.4)	GY-96410-17	25		GY-96403-17	500	
L/S 18	3/8 (9.5)	GY-96410-18	25		GY-96403-18	400	
L/S® High-performance Precision pump tubing							
L/S 15	3/16 (4.8)	GY-96410-15	25		GY-96403-15	400	
L/S 24	1/4 (6.4)	GY-96410-24	25		GY-96403-24	300	
L/S 35	3/8 (9.5)	GY-96410-35	25		GY-96403-35	300	
L/S 36	3/8 (9.5)	GY-96410-36	25		GY-96403-36	100	
I/P® Precision pump tubing							
I/P 26	1/4 (6.4)	GY-96410-26	25		GY-96403-26	200	
I/P 73	3/8 (9.5)	GY-96410-73	25		GY-96403-73	150	
I/P 82	1/2 (12.8)	GY-96410-82	25		GY-96403-82	100	
I/P® High-Performance Precision pump tubing							
I/P 70	3/8 (9.5)	GY-96510-70	10	—	—	—	—
I/P 88	1/2 (12.8)	GY-96510-88	10	—	—	—	—
I/P 89	5/8 (15.9)	GY-96510-89	10	—	—	—	—
B/T® PerfectPosition™ pump tubing							
B/T 87	3/8 (9.5)	GY-96510-87	10	—	—	—	—
B/T 91	3/4 (19.0)	GY-96510-91	10	—	—	—	—



Tubing size	Hose barb in. (mm)	Catalog number	Ft/pk	Price/pk
L/S® Precision pump tubing				
L/S 13	1/16 (1.6)	GY-96440-13	25	
L/S 14	1/16 (1.6)	GY-96440-14	25	
L/S 16	1/8 (3.2)	GY-96440-16	25	
L/S 25	3/16 (4.8)	GY-96440-25	25	
L/S 17	1/4 (6.4)	GY-96440-17	25	
L/S 18	3/8 (9.5)	GY-96440-18	25	
L/S High-Performance Precision pump tubing				
L/S 15	3/16 (4.8)	GY-96440-15	25	
L/S 24	1/4 (6.4)	GY-96440-24	25	
L/S 35	3/8 (9.5)	GY-96440-35	25	
L/S 36	3/8 (9.5)	GY-96440-36	25	
I/P® Precision pump tubing				
I/P 26	1/4 (6.4)	GY-96440-26	25	
I/P 73	3/8 (9.5)	GY-96440-73	25	
I/P 82	1/2 (12.8)	GY-96440-82	25	
B/T® PerfectPosition™ pump tubing				
B/T 87	3/8 (9.5)	GY-96445-87	10	
B/T 91	3/4 (19.0)	GY-96445-91	10	

Wait! There's More...

Masterflex® tubing formulation information, see pages 679-681.

Masterflex® Peroxide-Cured Silicone Pump Tubing

See samples better



Applications: Pharmaceutical, biotechnology, peristaltic pumping, and food applications.
Characteristics: Excellent heat resistance, low compression set, odorless, nontoxic. Translucent.
Certifications: USP Class VI, European Pharmacopoeia (EP 3.1.9), FDA 21 CFR 177.2600, exceeds 3A, REACH (non-DEHP) Compliant, RoHS Compliant, ADFC Compliant.
Temperature range: -58 to 446°F (-50 to 230°C)
Sterilize: Autoclave or ethylene oxide.



ANIMAL DERIVED
INGREDIENT FREE

Tubing size	Hose barb in. (mm)	Catalog number	Ft/pk	Price/pk
L/S® Precision pump tubing				
L/S 13	1/16 (1.6)	GY-96400-13	25	
L/S 14	1/16 (1.6)	GY-96400-14	25	
L/S 16	1/8 (3.2)	GY-96400-16	25	
L/S 25	3/16 (4.8)	GY-96400-25	25	
L/S 17	1/4 (6.4)	GY-96400-17	25	
L/S 18	3/8 (9.5)	GY-96400-18	25	
L/S® High-Performance Precision pump tubing				
L/S 15	3/16 (4.8)	GY-96400-15	25	
L/S 24	1/4 (6.4)	GY-96400-24	25	
L/S 35	3/8 (9.5)	GY-96400-35	25	
L/S 36	3/8 (9.5)	GY-96400-36	25	
I/P® Precision pump tubing				
I/P 26	1/4 (6.4)	GY-96400-26	25	
I/P 73	3/8 (9.5)	GY-96400-73	25	
I/P 82	1/2 (12.8)	GY-96400-82	25	
I/P® High-Performance Precision pump tubing				
I/P 70	3/8 (9.5)	GY-96400-70	10	
I/P 88	1/2 (12.8)	GY-96400-88	10	
I/P 89	5/8 (15.9)	GY-96400-89	10	
B/T® PerfectPosition™ pump tubing				
B/T 87	3/8 (9.5)	GY-96400-87	10	
B/T 91	3/4 (19.0)	GY-96400-91	10	