

### Cole-Parmer® Premium Luer Fittings

Ensure your process integrity with full material traceability

– Manufactured from 100% virgin raw materials in a Class 100,000 environment

These luer fittings are ideal for a variety of critical applications, including syringes or small tubing. Manufactured in compliance with ISO 594-1 standards, these fittings are inspected to be free from contaminants, and are fully traceable. All parts meet tight dimensional tolerances for increased reliability and process purity. Choose from a range of sizes, 1/16" to 1/4", with a variety of connection options. **Sold in pack of 25.**

Luer Fittings



Key	Description	Polypropylene		Polycarbonate		Nylon		PVDF	
		Catalog number	Price/pk	Catalog number	Price/pk	Catalog number	Price/pk	Catalog number	Price/pk
<b>Female luer x standard hose barb adapters</b>									
A	Female luer x 1/16" hose barb	<a href="#">GY-45508-00</a>		<a href="#">GY-45501-00</a>		<a href="#">GY-45502-00</a>		<a href="#">GY-45512-00</a>	
	Female luer x 3/32" hose barb	<a href="#">GY-45508-02</a>		<a href="#">GY-45501-02</a>		<a href="#">GY-45502-02</a>		<a href="#">GY-45512-02</a>	
	Female luer x 1/8" hose barb	<a href="#">GY-30800-08</a>		<a href="#">GY-45501-04</a>		<a href="#">GY-45502-04</a>		<a href="#">GY-45512-04</a>	
	Female luer x 5/32" hose barb	<a href="#">GY-45508-06</a>		<a href="#">GY-45501-06</a>		<a href="#">GY-45502-06</a>		<a href="#">GY-45512-06</a>	
	Female luer x 3/16" hose barb	<a href="#">GY-45508-08</a>		<a href="#">GY-45501-08</a>		<a href="#">GY-45502-08</a>		<a href="#">GY-45512-08</a>	
	Female luer x 1/4" hose barb	<a href="#">GY-30800-06</a>		<a href="#">GY-45501-20</a>		<a href="#">GY-45502-20</a>		—	
<b>Female luer x hose barb adapters (smooth bore)</b>									
B	Female luer x 1/16" hose barb	<a href="#">GY-45508-50</a>		<a href="#">GY-45501-50</a>		<a href="#">GY-45502-50</a>		—	
	Female luer x 3/32" hose barb	<a href="#">GY-45508-52</a>		<a href="#">GY-45501-52</a>		<a href="#">GY-45502-52</a>		—	
	Female luer x 1/8" hose barb	<a href="#">GY-45508-54</a>		<a href="#">GY-45501-54</a>		<a href="#">GY-45502-54</a>		—	
<b>Female luer x hose barb adapters (wide bore)</b>									
C	Female luer x 1/16" hose barb	<a href="#">GY-45508-01</a>		<a href="#">GY-45501-01</a>		<a href="#">GY-45502-01</a>		—	
	Female luer x 3/32" hose barb	<a href="#">GY-45508-03</a>		<a href="#">GY-45501-03</a>		<a href="#">GY-45502-03</a>		—	
	Female luer x 1/8" hose barb	<a href="#">GY-30800-02</a>		<a href="#">GY-45501-05</a>		<a href="#">GY-45502-05</a>		—	
<b>Female luer x hose barb adapters (low-profile barb)</b>									
D	Female luer x 1/16" hose barb	<a href="#">GY-45508-12</a>		<a href="#">GY-45501-12</a>		<a href="#">GY-45502-12</a>		—	
	Female luer x 3/32" hose barb	<a href="#">GY-45508-14</a>		—		<a href="#">GY-45502-14</a>		<a href="#">GY-45512-14</a>	
	Female luer x 1/8" hose barb	<a href="#">GY-45508-16</a>		<a href="#">GY-45501-16</a>		<a href="#">GY-45502-16</a>		<a href="#">GY-45512-16</a>	
	Female luer x 0.17" hose barb	<a href="#">GY-45508-18</a>		—		<a href="#">GY-45502-18</a>		<a href="#">GY-45512-18</a>	
<b>Female luer x thread adapters</b>									
E	Female luer x 10-32 special tapered thread	<a href="#">GY-45508-70</a>		<a href="#">GY-45501-70</a>		<a href="#">GY-45502-70</a>		—	
	Female luer x 10-32 UNF thread	<a href="#">GY-45508-60</a>		<a href="#">GY-45501-60</a>		<a href="#">GY-45502-60</a>		—	
	Female luer x 1/4-28 UNF thread	<a href="#">GY-45508-66</a>		<a href="#">GY-45501-66</a>		<a href="#">GY-45502-66</a>		<a href="#">GY-45512-66</a>	
	Female luer x 1/4-28 UNF bottom sealing thread	<a href="#">GY-45508-76</a>		—		<a href="#">GY-45505-76</a>		<a href="#">GY-45513-76</a>	
	Female luer x 1/8-27 NPT thread	<a href="#">GY-30800-00</a>		—		<a href="#">GY-45505-78</a>		<a href="#">GY-45513-78</a>	
F	Female luer x 1/8-27 NPT thread elbow	<a href="#">GY-45518-80</a>		—		<a href="#">GY-45505-80</a>		<a href="#">GY-45513-80</a>	
<b>Female luer bulkhead x hose barb adapters</b>									
G	Female luer x 1/4-28 UNF thread to 1/16" hose barb	<a href="#">GY-45508-30</a>		<a href="#">GY-45501-30</a>		<a href="#">GY-45502-30</a>		<a href="#">GY-45512-30</a>	
	Female luer x 1/4-28 UNF thread to 3/32" hose barb	<a href="#">GY-45508-32</a>		<a href="#">GY-45501-32</a>		<a href="#">GY-45502-32</a>		<a href="#">GY-45512-32</a>	
	Female luer x 1/4-28 UNF thread to 1/8" hose barb	<a href="#">GY-45508-34</a>		<a href="#">GY-45501-34</a>		<a href="#">GY-45502-34</a>		<a href="#">GY-45512-34</a>	
<b>Female adapters</b>									
H	Female luer x female luer	<a href="#">GY-45508-22</a>		<a href="#">GY-45501-22</a>		<a href="#">GY-45502-22</a>		<a href="#">GY-45512-22</a>	
I	Female tee	<a href="#">GY-45508-56</a>		<a href="#">GY-45501-56</a>		<a href="#">GY-45502-56</a>		—	
J	Female elbow	<a href="#">GY-45508-26</a>		<a href="#">GY-45501-26</a>		<a href="#">GY-45502-26</a>		<a href="#">GY-45512-26</a>	
K	Female luer x cap	<a href="#">GY-30800-12</a>		<a href="#">GY-45501-28</a>		<a href="#">GY-45502-28</a>		<a href="#">GY-45512-28</a>	
L	Female luer x cap with tether loop	<a href="#">GY-45508-64</a>		<a href="#">GY-45501-64</a>		<a href="#">GY-45502-64</a>		<a href="#">GY-45512-64</a>	
<b>Male x female adapters</b>									
M	Male luer x female luer coupler	<a href="#">GY-45508-82</a>		<a href="#">GY-45501-82</a>		<a href="#">GY-45502-82</a>		<a href="#">GY-45512-82</a>	
N	Male luer with lock ring x female luer coupler	<a href="#">GY-45508-80</a>		<a href="#">GY-45501-80</a>		<a href="#">GY-45502-80</a>		—	
O	Male luer x female luer elbow	<a href="#">GY-45508-84</a>		<a href="#">GY-45501-84</a>		<a href="#">GY-45502-84</a>		—	
P	Male x female x female tee	<a href="#">GY-45508-85</a>		<a href="#">GY-45501-85</a>		<a href="#">GY-45502-85</a>		—	
Q	Male x male x female tee	<a href="#">GY-45508-75</a>		<a href="#">GY-45501-75</a>		<a href="#">GY-45502-75</a>		—	